

**AquaUltra Clear Ultrasound Transmission Gel** has been designed for use in ultrasound physiotherapy and diagnostic applications. The gel guarantees perfect transmission of ultrasound waves.

- Water soluble, odourless, colourless, non-staining, non-irritant, easy to remove
- Contains no fatty substances, salt or formaldehyde
- Does not contain any active ingredients which may damage probes of echograph equipment

**INSTRUCTIONS FOR USE**

- The gel is ready for use
- Clean and degrease the patient's skin
- Apply the desired amount of gel
- Connect and use equipment and probes following manufacturer's instructions
- Any residue of the gel can be removed with a paper towel or gauze moistened with water
- Stop using the gel for hypersensitive patients, if skin reactions occur

**WARNINGS**

- External use only
- Product is non sterile
- Do not apply the gel on skin affected by abrasions, erythema or wounds
- Close the bottle after use
- May cause irritation or redness with allergic persons
- Contact with eyes may cause irritation
- Analyst with sensitive skin should wear protective gloves during use
- Do not use the gel after the "use by" date on the bottle

**INGREDIENTS:** aqua, carbomer, sodium hydroxide, EDTA, phenoxyethanol, ethylhexylglycerin (colourless)

**PACKAGING:** Carton contains 35 bottles of 260g













**DISPOSAL CONSIDERATIONS**

Product, waste: Chemicals must be disposed of in compliance with local and national regulations  
Packaging must be disposed of in compliance with local and national regulations

**STORAGE**

Temperature: +10°C - +30°C Relative humidity: 30% - 75%  
Keep container firmly closed, away from direct sunlight and sources of UV radiation

**EXPIRY:** 36 months from production date, use within 3 months of opening  
01.02.2021 V02

	Manufacturer		Temperature limit		Batch-Code		Attention
	User Manual		Polyethylene packaging		Use by date		Do not re-use
	Non sterile		Away from UV		Medical Device		

Manufacturer:  
**Ultrigel Medical Kft.**

✉ 1022 Budapest Aranka út 12. ☎ +36 1 278 3050  
📧 [info@ultrigel.hu](mailto:info@ultrigel.hu) 🌐 [www.ultrigel.hu](http://www.ultrigel.hu)

Distributor:  
**Ebrington (NI) Ltd**

✉ Maydown BT47 6UQ ☎ +44 (0)28 7134 5188  
📧 [info@ebrington-ni.com](mailto:info@ebrington-ni.com) 🌐 [www.eb-med.com](http://www.eb-med.com)