


product information

 Electric or Hydraulic height adjustment 45-96cm	 Maximum user weight (40st.) 250kg	 Gas strut assisted multi-position head and foot sections	 Back rest and leg rest adjustment +70° / -35°	 Easy to use handset (Electric models only)
 Retractable positioning castors	 Adjustable height levelling foot	 Fully upholstered head & foot sections	 Double layer density foam for support & comfort	 Luxury 65mm deep contoured upholstery
 Anti-bacterial/anti-microbial vinyl	 Paper roll holder (supplied as standard)	 Easy clean design and materials	 FIRA Independently tested by FIRA Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998	 West Yorkshire Materials Testing Service Anti-bacterial Vinyl Fire (Flammability) Testing BS 5852 : 2006 Clause 11 (composites)

accessories \* these accessories can only be fitted when ordered with the couch

 Sun-FECO/LEG1 Leg Supports	 Sun-FECO/PHB1 Phlebotomy Arm (each)	 Sun-FECO/FSW1* Foot Switch	 Sun-FECO/SUP1 Gynae Stirrups	 Sun-CHP1 Head Pillow
 Sun-FECO/CBH1* Breathing Hole	 Sun-FECO/CSR* Folding Side Rails (pair)	 Sun-FECO/IV2* IV Pole Attachment	 Sun-FECO/RAIL* Medi-Rail Attachment	 Sun-FECO/BB* Battery Backup
 Sun-FECO/HEO* Electrically Operated Head Section Upgrade	 Sun-FECO/FEO* Electrically Operated Foot Section Upgrade	 Sun-INST1* Installation Service & Packaging Recycling	 Sun-DISPOSAL1** Removal & Disposal of Old Couch	<b>** Removal from ground floor only - weight &amp; removal restrictions apply, please enquire</b>

dimensions




185cm(L) 65cm(W)



anti-bacterial vinyl upholstery colours - please specify when ordering


please see inside back cover for full specification of our vinyl range



The new Fusion Eco couch range provides most of the functionality you would expect for our standard Fusion range at a more affordable price to suit individual budgets. Fusion Eco is simply, strong, stable, aesthetically pleasing and easy to clean.



Fully upholstered head & foot sections with gas assisted adjustment



Electric or Hydraulic height adjustment



Increased maximum user weight to 250kg



**Optional**  
Electric Head & Foot Sections  
(electric models only)



Code: Sun-FECO3E  
Fusion Eco 3 Section Electric Couch (Shown with optional accessories)



3 year warranty on all moving components



Normally delivered within 7-10 working days (\*UK Only)



Lifetime warranty on frame



Code: Sun-FECO3H  
Fusion Eco 3 Section Hydraulic Couch



Hydraulic height adjustment