

product information

| | | | | |
|---|--|---|--|---|
| | | | | |
| Electric or Hydraulic height adjustment 45-98cm | Maximum user weight (35st.) 225kg | Gas assisted back (-20 - +80°) | Easy to use handset (Electric models only) | Rear retractable positioning castors |
| | | | | |
| 100mm Front braked castors | Cut away lithotomy section for easier access | 100mm Deep stainless steel operators tray | Fully adjustable leg supports | Fold down arms for patient comfort & security |
| | | | | |
| Double layer density foam for support & comfort | Anti-bacterial/anti-microbial vinyl | Easy clean design and materials | Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998 | Fire (Flammability) Testing Flammability BS 7176 Medium Hazard BS 5852: 2006 Clause 11 (composites) |

accessories

* these accessories can only be fitted when ordered with the couch

| | | | | |
|--|--|----------------------------------|--|---|
| | | | | |
| Sun-FGYN/FSW1* Foot Switch (For Gynae 1) | Sun-FGYN/CHP1 Head Pillow | Sun-FGYN/CBH1* Breathing Hole | Sun-FGYN/IV2* IV Pole Attachment | Sun-FGYN/SUP1 Stirrups (pair) |
| | | | | |
| Sun-FGYN/FS Detachable Foot Section Extension | Sun-FGYN/WB Wall Bracket for Foot Section | Sun-FGYN/BB* Battery Backup | Sun-INST1* Installation Service & Packaging Recycling | Sun-DISPOSAL1/** Removal & Disposal of Old Couch |

** Removal from ground floor only - weight and removal restrictions apply, please enquire

dimensions

137cm(L) 65cm(W)
187cm with foot extension

45-98cm(H)
-20 - +80°

Detachable Foot Section Extension
50cm(L) 65cm(W)

HEAVY ITEM

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range

| | | | | | | | | | | | | | | |
|-------|------|-------|----------|--------|----------|--------|-------|------|----------|----------|-----------|-------|------|-------|
| | | | | | | | | | | | | | | |
| black | grey | beige | primrose | walnut | red wine | salmon | lilac | navy | mid blue | sky blue | cool blue | green | mint | white |

The new Gynae couch range features an ergonomic design providing comfort for the patient and ease of access and patient positioning for the clinician. Fusion Gynae is simply, strong, stable, aesthetically pleasing and easy to clean.



Code: **Sun-FGYNE1**
Fusion Gynae1 - 2 section couch
Electric powered height, gas assisted back, no seat tilt with arm supports & leg supports

3 year warranty on all moving components

Lifetime warranty on frame



Normally delivered within 7-10 working days (*UK Only)

Code: **Sun-FGYNH1**
Fusion Gynae1 - 2 section couch
Hydraulic powered height, gas assisted back, no seat tilt with arm supports & leg supports



Hydraulic height adjustment

product information


| | | | | |
|--|--|--|--|--|
|  Electric height adjustment 45-98cm |  225kg Maximum user weight (35st.) |  Electric Powered Tilting Seat 0 to -15° (Trendelenburg) |  Fully powered head section (Sun-FGYNE3) |  Easy to use handset |
|  Rear retractable positioning castors |  100mm Front braked castors |  Cut away lithotomy section for easier access |  100mm Deep stainless steel operators tray |  Fully adjustable leg supports |
|  Fold down arms for patient comfort & security |  Double layer density foam for support & comfort |  Easy clean design and materials |  Independently tested by FIRA |  Anti-bacterial Vinyl |
| | | Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998 | | Fire (Flammability) Testing BS 5852: 2006 Clause 11 (composites) |

accessories

* these accessories can only be fitted when ordered with the couch

| | | | | |
|---|---|---|--|---|
|  Sun-FGYN/FSW2* Foot Switch (Gynae 2) |  Sun-FGYN/FSW3* Foot Switch (Gynae 3) |  Sun-FGYN/CHP1 Head Pillow |  Sun-FGYN/CBH1* Breathing Hole |  Sun-FGYN/IV2* IV Pole Attachment |
|  Sun-FGYN/SUP1 Stirrups (pair) |  Sun-FGYN/FS Detachable Foot Section Extension |  Sun-FGYN/WB Wall Bracket for Foot Section |  Sun-FGYN/BB* Battery Backup |  Sun-INST1* Installation Service & Packaging Recycling |
|  Sun-DISPOSAL1** Removal & Disposal of Old Couch | ** Removal from ground floor only - weight & removal restrictions apply, please enquire | | | |

dimensions



137cm(L) 65cm(W)
187cm with foot extension

50cm(L) 65cm(W)

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range

| | | | | | | | | | | | | | | |
|---|---|---|---|---|---|---|---|---|---|---|--|---|---|---|
|  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| black | grey | beige | primrose | walnut | red wine | salmon | lavender | navy | mid blue | sky blue | cool blue | green | mint | white |

The new Gynae couch range features an ergonomic design providing comfort for the patient and ease of access and patient positioning for the clinician. Fusion Gynae is simply, strong, stable, aesthetically pleasing and easy to clean.



Code: Sun-FGYNE2
Fusion Gynae2 - 2 section couch
Electric powered height, gas assisted back, powered tilting seat with arm supports & leg supports



3 year warranty on all moving components



Lifetime warranty on frame



Normally delivered within 7-10 working days (*UK Only)



Code: Sun-FGYNE3
Fusion Gynae3 - 2 section couch
Electric powered height, back and tilting seat with arm supports & leg supports

