

product information


Maximum user weight (17st.)	Lightweight design only 8kgs	Multi-position arm rest, made from a hygienic PU moulding	Arm rest can be positioned either side (2 brackets supplied)	Ergonomically designed moulded seat & back (Sun-PCHA only)
Anti-bacterial/anti-microbial vinyl (Sun-PCHA/VYL only)	40mm Deep foam for support and comfort (Sun-PCHA/VYL only)	Easy clean design and materials	1 year warranty	Supplied fully assembled
Independently tested	Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998	Fire (Flammability) Testing Flammability BS 7176 Medium Hazard BS 5852 : 2006 Clause 11 (composites)	British made	CE marked

accessories



Sun-PCHA/Arm
Additional Arm Rest

dimensions



83cm(H) 70cm(W) 56cm(D)
Including arm brackets (excluding arm)
Seat Height 43.5cm(H)

Arm Rest
35cm(D)
49cm(H)

Maximum user weight 110kg (17st.)

12cm(D)
35cm(H)

moulded seat & back colours - please specify when ordering



anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



Durable, stable and ergonomically designed, phlebotomy chairs are easy to use and adjust and can be used left, right or dual sided with an additional arm, whilst also being easy to clean.



Code: **Sun-PCHA**
Phlebotomy chair
Moulded seat and back, one arm rest



CODE: **Sun-PCHA/VYL**
Phlebotomy chair
Vinyl upholstery, one arm rest

