

## Fusion Bariatric Couches

Our NEW Fusion Bariatric Couch extends our already popular couch range. With on going requirements for examination couches to have the capability to support increasingly heavy patients we have designed our 75cm wide bariatric couch to have a safe working load/lifting capacity of 320kg (50st.) Provides enhanced user functionality and easier patient positioning through its fully electric powered height and head section. Our Fusion Bariatric Couch is simply, strong, stable, aesthetically pleasing and easy to clean.

- Maximum user weight 320kg (50st.)
- Available in 2 or 3 section
- 75cm wide for easier patient positioning (100cm wide also available)
- Electric actuators for main up/down movement and head section
- Easy to use handset
- 4 Large positioning castors with central locking system
- Fully upholstered head & foot sections
- Double layer density foam for support & comfort
- Luxury 65mm deep contoured upholstery
- 15 Anti-bacterial/anti-microbial vinyl colours available
- Anti-bacterial Vinyl - Fire (Flammability) Testing  
BS 7176 Medium Hazard BS 5852 : 2006 Clause 11
- Paper roll holder (supplied as standard)
- Silver powder coated frame
- Easy clean design and materials
- 3 year warranty on all moving components



Sun-FBVHC2E/75




Sun-FBVHC3E/75




<b>Sun-FBVHC2E/75</b> Fusion Bariatric 2 Section Variable Height Electric Couch	<b>£1250.00</b>
<b>Sun-FBVHC3E/75</b> Fusion Bariatric 3 Section Variable Height Electric Couch	<b>£1450.00</b>

### dimensions




### accessories

(\*These accessories can only be fitted when ordered with couch)




**Sun-FBVHC/CBH1\***  
Breathing Hole  
£35.00



**Sun-FBVHC/CHP1**  
Head Pillow  
£75.00




**Sun-FBVHC/FSW2\***  
Foot Switch for 2 Section  
£215.00



**Sun-FBVHC/FSW3\***  
Foot Switch for 3 Section  
£315.00



**Sun-FBVHC/CSR\***  
Folding Side Rails (pair)  
£175.00



**Sun-FBVHC/BB\***  
Battery Backup  
£149.00



**SUN-DISPOSAL1\***  
Removal and Disposal  
of Old Couch  
£95.00

### anti-bacterial vinyl upholstery colours - please specify when ordering

