


product information

Electric or Hydraulic height adjustment	Maximum user weight (40st.)	50mm Retractable positioning castors	Gas assisted head and split foot sections	Electric head section adjustment (Sun-FPODE2)
Adjustable leg length extension upto 15cm	Tilting Seat for Trendelenburg (Sun-FPODE1 & Sun-FPODE2)	Stay level swing back padded arm rests for easy patient access	Adjustable height levelling foot	Hoist clearance 13cm
Fully upholstered head & foot sections	Double layer density foam for support & comfort	Luxury 65mm deep contoured upholstery	Anti-bacterial/anti-microbial vinyl	Paper roll holder (supplied as standard)
Easy clean design and materials	Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998	Fire (Flammability) Testing Flammability BS 7176 Medium Hazard BS 5852 : 2006 Clause 11 (composites)	Anti-bacterial Vinyl	

accessories * these accessories can only be fitted when ordered with the couch


Sun-CHP1 Head Pillow	Sun-FPOD/IV2* IV Pole Attachment	Sun-FPOD/FSW1* Foot Switch (FPODE1)	Sun-FPOD/FSW2* Foot Switch (FPODE2)	Sun-FPOD/FSW2* Foot Switch (FPODE2)	installation
	** Removal from ground floor only - weight & removal restrictions apply, please enquire				
Sun-DISPOSAL1/** Removal & Disposal of Old Couch					

dimensions



anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



The new Fusion Podiatry Chair range has been designed using our extensive knowledge to provide a trusted, ergonomic design providing comfort for the patient, ease of access via the split leg design, Trendelenburg movement as standard on electric models and comfortable patient positioning for the clinician.



Code: **Sun-FPODE1**
Fusion Podiatry Chair
Electric Height Adjustment
Gas assisted head and split adjustable length foot sections
(Powered Trendelenburg -20° Movement)

Code: **Sun-FPODE2** (shown above)
Fusion Podiatry Chair
Electric Height Adjustment
Electric head and split adjustable length foot sections
(Powered Trendelenburg -20° Movement)



3 year warranty on all moving components



Lifetime warranty on frame



Normally delivered within 7-10 working days (*UK Only)

Code: **Sun-FPODH1**
Fusion Podiatry Chair
Hydraulic Height Adjustment
Gas assisted head and split adjustable length foot sections
(No Trendelenburg Movement)



Hydraulic height adjustment