

product information

 Maximum user weight (26st.)	 Multi-position self locking backrest	 Single layer density foam for support & comfort	 Luxury 65mm deep contoured upholstery	 Anti-bacterial/anti-microbial vinyl
 Paper roll holder (supplied as standard)	 Durable silver epoxy coated frame	 Stylish 60cm diameter legs for stability	 4 non marking leg protectors	 Special order upholstery colours available
 Easy clean design and materials	 Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998	 Fire (Flammability) Testing Flammability BS 7176 Medium Hazard BS 5852 : 2006 Clause 11 (composites)		

accessories * these accessories can only be fitted when ordered with the couch

 Sun-CBH1* Breathing Hole	 Sun-CHP1 Head Pillow	 Sun-LEG1* Leg Supports (pair)	 Sun-INST1* Installation Service & Packaging Recycling	 Sun-DISPOSAL1*/** Removal & Disposal of Old Couch
 Sun-MPT51 Vista 15 Under Couch Trolley	 Sun-MPT52 Vista 15 Under Couch Trolley	 Sun-MPT53 Vista 15 Under Couch Trolley	 Vista 15 Under Couch Trolleys (see pages 56-57)	 Couch Steps (see pages 42-43)

** Removal from ground floor only - weight and removal restrictions apply, please enquire

dimensions



anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



Simple, strong and durable, Colenso Couches are designed for everyday general healthcare, providing comfort with contemporary design.



Code: **Sun-MPT51**
Vista 15 under couch trolley
4 single depth trays



Code: **Sun-MPT52**
Vista 15 under couch trolley
2 single/1 double depth trays



Code: **Sun-MPT53**
Vista 15 under couch trolley
2 double depth trays