


product information

 Powered tilting seat Sun-FTRTE2/3/5/6 only	 -20° TRENDELEBURG Tilting Seat for Trendelenburg Sun-FTRTE2/3/5/6 only	 Head adjustment +80° / -20°	 Adjustable foot section +20° / -90°	 2 Adjustable/removable padded arm supports
 250kg Maximum user weight (40st.)	 Easy to use handset	 50mm Retractable positioning castors	 Adjustable height levelling foot	 Dual sided foot pedals (Hydraulic models only)
 Hoist clearance 13cm	 Fully upholstered head & foot sections	 Double layer density foam for support & comfort	 Luxury 65mm deep contoured upholstery	 Anti-bacterial/ anti-microbial vinyl
 Paper roll holder (supplied as standard)	 Easy clean design and materials	 Independently tested by FIRA Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998	 West Yorkshire Materials Testing Service Anti-bacterial Vinyl Fire (Flammability) Testing Flammability BS 7176 Medium Hazard BS 5852 : 2006 Clause 11 (composites)	

accessories * these accessories can only be fitted when ordered with the couch

 Sun-FTRT/PHB1 Phlebotomy Arm (pair)	 Sun-CHP1 Head Pillow	 Sun-FTRT/CBH1* Breathing Hole - in foot section (Sun-FTRTE1,2,3 & Sun-FTRTH1 only)	 Sun-FTRT/FSW1* Foot Switch (Sun-FTRTE1 & 4)	 Sun-FTRT/FSW2* Foot Switch 2 Section (Sun-FTRTE2 & 5)
 Sun-FTRT/FSW3* Foot Switch 3 Section (Sun-FTRTE3 & 6)	 Sun-FTRT/IV2* IV Pole Attachment	 Sun-INST1* Installation Service & Packaging Recycling	 Sun-DISPOSAL1** Removal & Disposal of Old Couch	** Removal from ground floor only - weight & removal restrictions apply, please enquire

dimensions



170cm(L) 65cm(W)

**HEAVY
ITEM**

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



The new Fusion Treatment Chair range has been designed using our extensive knowledge to provide a trusted, ergonomic design providing comfort for the patient, ease of access via the split leg design and comfortable patient positioning for the clinician.



Powered tilting seat



Adjustable foot section



Recovery position



Code: **Sun-FTRTE3**
Fusion Treatment Chair - Electric Height Single Foot Section
Powered head section and powered tilting seat

Electric Height Single Foot Section

Code: **Sun-FTRTE1** - Gas assisted head section and fixed seat
Code: **Sun-FTRTE2** - Gas assisted head section and powered tilting seat
Code: **Sun-FTRTE3** - Powered head section and powered tilting seat

Electric Height Split Foot Section

Code: **Sun-FTRTE4** - Gas assisted head section and fixed seat
Code: **Sun-FTRTE5** - Gas assisted head section and powered tilting seat
Code: **Sun-FTRTE6** - Powered head section and powered tilting seat

Hydraulic Height Adjustment with Fixed Seat

Code: **Sun-FTRTH1** - Gas assisted head section and single foot section
Code: **Sun-FTRTH4** - Gas assisted head section and split foot section



3 year warranty on all
moving components



Lifetime warranty
on frame



Normally delivered within
7-10 working days (*UK Only)



Code: **Sun-FTRTH4**
Fusion Treatment Chair - Hydraulic Height Adjustment Fixed Seat
Gas assisted head section and split foot section

