

VACUAIDE 7325 PORTABLE SUCTION UNIT

For effective
airways clearance

Very
lightweight
at only 1.7kg
including the
battery

Li-ion battery
provides a fast
charging time
and a long
battery life



The 7325P is
supplied with its
own carry bag
designed to allow
discreet portable use
of the device



For convenience
and hygiene an
integrated tube holder
allows the user to
place the suction tube
after each suction
procedure



The VacuAide 7325 Portable Suction Unit has been designed to perform effective airways suctioning in a comfortable, safe and hygienic way. Featuring a lightweight Li-ion battery with a quick charging time and long battery life for improved performance. The VacuAide 7325 is suitable for both a clinical and homecare setting.

www.drivedevilbiss.co.uk

drive
DeVilbiss[®]
HEALTHCARE



Features:

- » Features an internal Li-ion battery (7325P) with a fast charge time of 3 hours, a runtime of 45 minutes and a 1 year warranty on the battery
- » Very lightweight at 1.7kg including the battery for easy transportation
- » A suction tube holder is integrated into the 7325 for convenience and hygiene
- » The vacuum gauge is shock protected and can be adjusted from 50 - 550mmHg
- » Easy to operate and maintenance free
- » Can be operated discreetly inside the personal carry case
- » The personal carry case includes a viewing panel so that you can visually check the secretion level at any time. It is padded and is designed as a shoulder bag, with internal compartments for accessories
- » Complete with 100-240 V power supply and 12 V DC car adaptor
- » Collection canister is easy to install and remove
- » The VacuAide 7325 comes complete with an 800ml collection canister
- » A reusable autoclavable 1200ml canister and the Flovac canister and disposable liner system are available as accessories
- » LED control lights are displayed on the top of the device for quickly checking the units charge status
- » High suction power – up to 27 l/min free flow
- » To prevent liquids entering the suction unit there are two safety features; a hydrophobic filter and a mechanical overflow protection built into the canister lid
- » Shares many of the consumables with the existing VacuAide suction devices, including collection canisters, filters and patient tubing
- » A mains powered version is also available (D model). Please note this doesn't include the internal battery, 12 V car adapter or carry case

Product Codes:

Code	Description
7325P-U	VacuAide 7325 Portable Suction Unit (with internal rechargeable battery)
7325D-U	VacuAide 7325 Mains Powered Suction Unit (without internal battery, 12V car adapter or carry case)

Specifications:

Height: 17.5cm	Width: 17.8cm
Depth: 18cm	Operating Temperature Range: 0°C - 40°C
Product Weight 7325P-U: 1.7kg	Product Weight 7325D-U: 1.5kg
Air Flow: 27 l/min	Vacuum Range: 50 – 550mmHg
Charging time (for 7325P): Up to 3 hours	Battery running time (for 7325P): Up to 45 mins
Device Warranty: 2 Years	Battery Warranty: 1 Year
Electrical Requirements: 100-240V AC, 50/60 Hz; 12V DC	



If you would like to view this product, please contact your local stockist.



LL507

www.drivedevilbiss.co.uk
0845 0600 333