






product information

| | | | | |
|---|--|---|---|--|
|  Electric or Hydraulic height adjustment |  Maximum user weight (40st.) |  Gas assisted head & foot sections |  Electric head section & foot sections adjustment option (Sun-FPHBE2 & 3) |  Tilting Seat for Trendelenburg (Sun-FPHBE3) |
|  Two fully adjustable phlebotomy arms as standard |  50mm Retractable positioning castors |  Adjustable height levelling foot |  Dual sided foot pedals (Hydraulic models only) |  Hoist clearance 13cm |
|  Fully upholstered head & foot sections |  Double layer density foam for support & comfort |  Luxury 65mm deep contoured upholstery |  Anti-bacterial/anti-microbial vinyl |  Paper roll holder (supplied as standard) |
|  Easy clean design and materials |  Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998 |  Fire (Flammability) Testing Flammability BS 7176 Medium Hazard BS 5852 : 2006 Clause 11 (composites) | | |


accessories

* these accessories can only be fitted when ordered with the couch

| | | | | |
|--|--|--|--|---|
|  Sun-CHP1 Head Pillow |  Sun-FPHB/FSW1* Foot Switch (Sun-FPHBE1) |  Sun-FPHB/FSW3* Foot Switch 3 Section (Sun-FPHBE2) |  Sun-FPHB/FSW4* Foot Switch 4 Section (Sun-FPHBE3) |  Sun-FPHB/IV2* IV Pole Attachment |
|  Sun-INST1* Installation Service & Packaging Recycling |  Sun-DISPOSAL1** Removal & Disposal of Old Couch | | | |

** Removal from ground floor only - weight & removal restrictions apply, please enquire

dimensions




170cm(L) 65cm(W)

HEAVY ITEM

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



The new Fusion Phlebotomy Chair range has been designed using our extensive knowledge to provide a trusted, ergonomic design providing comfort for the patient and comfortable patient positioning for the clinician.



Easy removal of the phlebotomy arms allows the chair to be used as a standard drop end couch



In emergency situations the chair goes head down and feet up for patient recovery



Code: **Sun-FPHBE3**
Fusion Phlebotomy Chair - Electric Height Adjustment
Electric back and foot sections and tilting seat

Electric Height Adjustment

Code: **Sun-FPHBE1** - Electric Height Adjustment gas assisted head and foot sections
Code: **Sun-FPHBE2** - Electric Height Adjustment electric back and foot sections
Code: **Sun-FPHBE3** - Electric back and foot sections and tilting seat



3 year warranty on all moving components



Lifetime warranty on frame



Normally delivered within 7-10 working days (*UK Only)



Code: **Sun-FPHBH1**
Fusion Phlebotomy Chair - Hydraulic Height Adjustment
Gas assisted head and foot sections



The mark of responsible forestry