


# Technical Data Sheet

## Button Tourniquet



<b>Product</b>	Single-use button tourniquet
<b>Branding</b>	Emerald Medical - Button Tourniquet
<b>Product Ref</b>	P272
<b>Primary Use</b>	Single-use, quick-release disposable tourniquet for restricting blood flow during blood sampling procedures
<b>Material</b>	Perforated thermoplastic elastomer band, 36cm x 2.5cm with rigid plastic 'mushroom' locking button
<b>Packing &amp; Labelling</b>	<p>Box contains 100 pieces. Marked with product name and description, product reference code, brand logo, GTIN data matrix barcode, CE marking, MD, single-use, Latex-free and storage requirement pictograms</p> <p>Outer carton contains 10 boxes. (Total product quantity 1000 pieces) Marked with product name and description, brand logo, product reference code, batch and lot identifiers, CE marking, MD, single-use, Latex-free and storage requirement pictograms</p>
<b>Storage Conditions</b>	Store in a cool dry place away from heat and direct sunlight
<b>Shelf Life</b>	3 years (if stored as directed above)
<b>Additional Information</b>	<p>Latex-free</p> <p>Single-use only</p>
<b>Classification</b>	Class I Medical Device Regulations 2017/745