


Technical Data Sheet

Paediatric Button Tourniquet



| | |
|--------------------------------|--|
| Product | Single-use paediatric button tourniquet |
| Branding | Emerald Medical - Paediatric Button Tourniquet |
| Product Ref | T069 |
| Primary Use | Single-use, quick-release disposable tourniquet for restricting blood flow during blood sampling procedures |
| Material | Perforated thermoplastic elastomer band, 30cm x 2cm with rigid plastic 'mushroom' locking button |
| Packing & Labelling | <p>Box contains 100 pieces. Marked with product name and description, product reference code, brand logo, GTIN data matrix barcode, CE marking, MD, single-use, Latex-free and storage requirement pictograms</p> <p>Outer carton contains 10 boxes. (Total product quantity 1000 pieces) Marked with product name and description, brand logo, product reference code, batch and lot identifiers, CE marking, MD, single-use, Latex-free and storage requirement pictograms</p> |
| Storage Conditions | Store in a cool dry place away from heat and direct sunlight |
| Shelf Life | 3 years (if stored as directed above) |
| Additional Information | <p>Latex-free</p> <p>Single-use only</p> |
| Classification | Class I Medical Device Regulations 2017/745 |