


Technical Data Sheet

Fishbone Tourniquet



Product	Single-use pull-through tourniquet (Fishbone)
Branding	Emerald Medical - Fishbone Tourniquet
Product Ref	T075
Primary Use	Single-use, pull-through disposable tourniquet for restricting blood flow during blood sampling procedures
Material	Thermoplastic profile-cut elastomer band 45cm x 2.5cm with pull-through apertures
Packing & Labelling	<p>Box contains 100 pieces. Marked with product name and description, product reference code, brand logo, GTIN data matrix barcode, CE marking, MD, single-use, Latex-free and storage requirement pictograms</p> <p>Outer carton contains 20 boxes. (Total product quantity 2000 pieces) Marked with product name and description, brand logo, product reference code, batch and lot identifiers, CE marking, MD, single-use, Latex-free and storage requirement pictograms</p>
Storage Conditions	Store in a cool dry place away from heat and direct sunlight
Shelf Life	3 years (if stored as directed above)
Additional Information	<p>Latex-free</p> <p>Single-use only</p>
Classification	Class I Medical Device Regulations 2017/745