



# STABIL

BUILT LIKE A ROCK



## PRO 2-SECTION TREATMENT TABLE

### ABOUT

The STABIL® name, synonymous with strength and quality is becoming the new standard for treatment tables worldwide.

The PRO series features angled legs and 4-point lift to provide ultimate stability while a closed base system and covered hinges improve hygiene. Ideal for busy therapy clinics where aesthetic appeal is key.

### FEATURES

- 2-Section Upholstery
- Luxury 50mm Foam
- Retractable Wheels
- Electric Hand Controller or Hydraulic Foot Pump
- Retractable Padded Side Railings
- Safe Working Load: 250KG



CAPACITY  
250kg



LENGTH  
195cm



WIDTH  
65cm

PRODUCT	SKU	ACTUATOR	LENGTH	WIDTH	HEIGHT	MAX WEIGHT	STRUCTURAL WARRANTY
2-SECTION / WHITE FRAME / BLUE UPHOLSTERY	1002H-65WBL 1002U-65WN	ELECTRIC / HYDRAULIC	195cm	65cm	43cm to 97cm	250kg	10 YEARS
2-SECTION TREATMENT TABLE / BLACK / 360 DEGREE OPERATION	1002U-65BB-ARB	ELECTRIC	195cm	65cm	43cm to 97cm	250kg	10 YEARS
2-SECTION - WHITE FRAME / BLUE UPHOLSTERY - 360 DEGREE OPERATION	1002U-65 WN-ARB	ELECTRIC	195cm	65cm	43cm to 97cm	250kg	10 YEARS
2-SECTION / BLACK FRAME / BLACK UPHOLSTERY	1002U-65BB	ELECTRIC	195cm	65cm	43cm to 97cm	250kg	10 YEARS

**OUR WARRANTY:** STABIL® warrants the Pro range to conform to specification and be free from defects in material and workmanship for a period of 10 years from date of purchase. Warranty covers frame and moving components (5 years) in normal service and excludes wear and tear to fabric & cables. Details and specification subject to change without notice E&OE.