

STABIL

BUILT LIKE A ROCK



KOMFORT 2-SECTION TREATMENT TABLE

ABOUT

The STABIL® name, synonymous with strength and quality is becoming the new standard for treatment tables worldwide.

Komfort™ combines strength and durability with comfortable, hard-wearing upholstery to enhance your treatment area for less.

FEATURES

- 2-Section Upholstery
- Firm 40mm Foam
- Retractable Wheels
- Electric Hand Controller or Hydraulic Foot Pump
- Optional Retractable Padded Side Railings
- Optional Gyne Stirrups



PRODUCT	SKU	ACTUATOR	LENGTH	WIDTH	HEIGHT	MAX WEIGHT	WARRENTY
2-SECTION / WHITE FRAME / NAVY UPHOLSTERY	SP132E SP132H	ELECTRIC / HYDRAULIC	191cm	63cm	46cm to 105cm	225kg	5 YEARS
2-SECTION / WHITE FRAME / BLACK UPHOLSTERY	SP232E SP232H	ELECTRIC / HYDRAULIC	191cm	63cm	46cm to 105cm	225kg	5 YEARS
2-SECTION / WHITE FRAME / MID BLUE UPHOLSTERY	SP332E SP332H	ELECTRIC / HYDRAULIC	191cm	63cm	46cm to 105cm	225kg	5 YEARS

OUR WARRENTY: STABIL® warrants the Komfort range to conform to specification and be free from defects in material and workmanship for a period of 5 years from date of purchase. Warranty covers frame and moving components in normal service and excludes wear and tear to fabric & cables. Details and specification subject to change without notice E&OE.

BUILT LIKE A ROCK™

Issue 10/23