




# STABIL

BUILT LIKE A ROCK



## KOMFORT 3-SECTION TREATMENT TABLE

### ABOUT

The STABIL® name, synonymous with strength and quality is becoming the new standard for treatment tables worldwide.

Komfort™ combines strength and durability with comfortable, hard-wearing upholstery to enhance your treatment area for less.

### FEATURES

- 3-Section Upholstery
- Firm 40mm Foam
- Retractable Wheels
- Electric Hand Controller or Hydraulic Foot Pump
- Optional Retractable Padded Side Railings
- Optional Gyne Stirrups



CAPACITY  
**225kg**



LENGTH  
**191cm**



WIDTH  
**63cm**

PRODUCT	SKU	ACTUATOR	LENGTH	WIDTH	HEIGHT	MAX WEIGHT	WARRENTY
3-SECTION / WHITE FRAME / NAVY UPHOLSTERY	SP133E SP133H	ELECTRIC / HYDRAULIC	191cm	63cm	46cm to 105cm	225kg	5 YEARS
3-SECTION / WHITE FRAME / BLACK UPHOLSTERY	SP233E SP233H	ELECTRIC / HYDRAULIC	191cm	63cm	46cm to 105cm	225kg	5 YEARS
3-SECTION / WHITE FRAME / MID BLUE UPHOLSTERY	SP333E SP333H	ELECTRIC / HYDRAULIC	191cm	63cm	46cm to 105cm	225kg	5 YEARS

**OUR WARRENTY:** STABIL® warrants the Komfort range to conform to specification and be free from defects in material and workmanship for a period of 5 years from date of purchase. Warranty covers frame and moving components in normal service and excludes wear and tear to fabric & cables. Details and specification subject to change without notice E&OE.