



**Ebrington**  
MEDICAL



## August/September Offers

### Stabil Pro 3-Section Treatment Table - Black Frame & Upholstery

Code: T958-ELE


Sturdy treatment table with angled legs ensure ultimate stability.

- Ultra Low Entry Height: 43cm
- Max Height: 97cm
- Weight Capacity: 250kgs
- Length: 195cm
- Premium quality upholstery
- Closed base system to hide all the electrical components

Was €1,492.56

Now €1,374.25

Offers in this brochure valid until 30/09/24



*We'll look after you!*

 [www.eb-med.com](http://www.eb-med.com)

T: 028 7134 5188

## seca CT331 PC Based Interpretive ECG Machine With Bluetooth

The innovative seca CT331 offers the best of both: Medical Bluetooth connectivity and USB cabled option.

Enjoy the flexibility of use in large clinic rooms with a wireless connection, but with a USB cabled option to plug in when needed. Cable-free operation also grants a higher freedom of movement and reduces trip hazards.

- Wireless and lightweight
- Transfer data via Bluetooth or USB cable
- Rechargeable inbuilt batteries for up to 12 hours of runtime (charge via USB port)
- Brilliant recording quality of all 12 leads with a sampling rate of up to 32 kHz
- LEDs for ECG quality and individual lead-on electrode check
- Clear surface for quick and easy disinfection
- EMR integration




Now **€2,300.38**

Code	Product	Price
T629	seca CT331 PC ECG with Bluetooth	€2,388.81



## seca CT8000P-2 Interpretive A4 ECG

- Robust and reliable: ideal for busy clinics
- Wi-Fi connectivity and seca link software included as standard – PDF ECGs automatically download to your PC following acquisition
- 8" colour high-resolution screen with anatomical electrode placement (hook-up adviser function)
- Compatible with all leading EMR systems
- Capture and store 10 minutes of uninterrupted 12 lead data to document intermittent arrhythmias
- Lead reversal: Warnings and misplaced electrodes are immediately displayed before ECG acquisition



Now **€3,579.03**

Code	Product	Price
M890	seca CT8000P-2 ECG	€3,848.03

## Welch Allyn CP150 ECG Machine with Interpretation

The extremely easy to use CP150 ECG is perfect for any fast-paced environment, where speed, simplicity and accuracy are all a priority.

### Features:

- 3 lead ECG screen preview
- Large, 7" colour touch-screen
- On-screen lead guidance
- A4 thermal printer for hard copy readings
- Easy to clean - reduce cross contamination
- Easily share or import results into your EMR
- Instant power up - ECGs at a moment's notice



Now **€2,802.28**

Code	Product	Price
T535	Welch Allyn CP150 ECG	€2,950.46

## Kendall Q Trace Gold Electrodes x 500

- Cardinal Health™ Kendall™ 5500 Diagnostic Tab Electrodes
- High-tack adhesive tab
- Designed for 12-lead ECG monitoring
- 100 pieces per pouch, 5 pouches



Now €29.28

Code	Product	Price
M325	Kendall Q Trace Electrodes	€34.66

## Skintact ECG Electrodes x 500

- Skintact EASITABS RT74
- Designed for resting 12 lead diagnostic ECGs
- 100 pieces per pouch, 5 pouches per box
- Excellent signal conduction



Now €29.28

Code	Product	Price
K382	Skintact ECG Electrodes	€34.66

## Ambu Blue Sensor Event Recording Electrodes x 25

- Ambu® Blue Sensor VLC
- Suitable for long-term or ambulatory monitoring
- Highly conductive wet gel
- Stud connector
- Breathable micro-pore material
- Pouch of 25



Now €22.71

Code	Product	Price
T281	Ambu Event Electrodes	€27.19

## Skintact ECG Monitoring Electrodes x 30

- Skintact T-V001
- For stress testing and Holter monitoring
- Reinforced adhesive
- PE foam backing
- Skin-friendly
- Stud connector
- Pouch of 30



Now €20.32

Code	Product	Price
K381	Skintact Monitoring Electrodes	€27.19

## ECG Skin Scratch Pads x 100

- Designed for skin preparation before electrode placement
- Abrasive to reduce skin impedance
- Adhesive on one side, to use with one finger
- 10 sheets of 10 pads



Now €22.71

Code	Product	Price
K479	Skin Scratch Pads	€30.83

## Snap Clip Connectors x 10

- Recommended for maximum connectivity to tab electrodes
- Resuable clip that connects 4 mm banana pins
- Ensures secure connection and are easy to use
- Pack of 10




Now €28.68

Code	Product	Price
K478	Snap Clip Connectors	€35.79

## Skin Prep Wipes x 100

- Swab saturated with 70% isopropyl alcohol v/v
- Sterile, individually wrapped
- 100 sachets per box



Now €1.55

Code	Product	Price
M166	Skin Prep Wipes	€2.37

## Disposable Prep Razors x 10

- Designed for the removal of body hair prior to patient surgery or monitoring
- Pack of 10 razors



Now €2.15

Code	Product	Price
M329	Disposable Razors	€2.63



## seca CT500-2 Height Adjustable ECG Trolley

The CT500-2 is a height adjustable, tip-proof ECG Cart, that provides easy transport for nearly all ECG types.

- Foot pedal to raise or lower the table effortlessly
- Storage bin with dividers for accessories
- Easy to clean and disinfect
- Braked castors for more stability and safety

Now **€543.73**


Code	Product	Price
M892	seca CT500-2 ECG Trolley	€594.53



## Fusion Treatment Chair with Single Foot Section - Electric Height and Back with Tilting Seat

Designed for patient comfort and ease of access for the clinician.

- Powered head section and powered tilting seat
- Two adjustable/removable padded arm supports as standard
- Maximum user weight 250kg
- Height range: 46-100cm
- Head adjustment +80°/-20°
- Foot adjustment +20°/-90°
- Choose from a range of colours



Now **€1,971.75**

Code	Product	Price
T806	Fusion Treatment Chair	€2,097.23



# LEC FRIDGE CLEARANCE SALE

**Grab a bargain.....once its gone, its gone!**

Manufacture Clearance sale - Lec Medical are looking to clear their excess stock of superseded pharmacy fridge models.

- Perfectly fit for use, never used and in original packaging
- As with any fridge, 2 year warranty starts from date of purchase
- Range of countertop and undercounter fridges
- Glass door or solid door options

Available Models Include:

- PPSR151UK            151L with solid door
- PPSR158UK-DLK    158L with solid door and digital lock
- PPSR47UK            47L with solid door
- PPGR47UK            47L with glass door
- PEGR107-LHH        107L with glass door and left hinge

**Contact us today for a quotation!**



# Coffee Break Challenge.....Matchstick Equations!!

Move only 1 matchstick to fix the equations



1)  $1 + 2 = 8$

Answer: .....

2)  $3 + 3 = 5$

Answer: .....

3)  $7 + 7 = 0$

Answer: .....

4)  $8 - 8 = 8$

Answer: .....

5)  $6 + 4 = 4$

Answer: .....

6)  $3 - 7 = 4$


Answer: .....

7)  $6 + 7 = 9$

Answer: .....

8)  $4 - 5 = 8$

Answer: .....

 **CLAIM YOUR PRIZE!**  
Send a copy of your answers to [customerservice@ebrington-ni.com](mailto:customerservice@ebrington-ni.com) with a note of which Practice you're from. Submit your correct entry by 16/08/24 and receive a small prize with your next order.